


Gold Nanocubes

DESCRIPTION

Biotyscience provides Gold Nanocubes of various sizes to meet the needs of different users for different purposes. Biotyscience's Gold Nanocubes products use deionized water as the solvent, and the CTAB content of the cosolvent does not exceed 0.01 M.



The product has a high degree of stability in water and can be evenly dispersed for several months without sedimentation. Each batch of products is synthesized through precise steps and characterized by independent visible/infrared extinction spectroscopy. Different batches of products have excellent repeatability.

Storage

Store product away from direct sunlight at 2-8 ° C.

Do NOT freeze. Freezing causes irreversible aggregation of the gold nanoparticles.

When stored as specified the product is stable for at least six months

PRODUCT INFORMATION

Type	Gold Nanocubes
Diameter	50 nm, 90 nm, 170 nm
Surface	CTAB
OD	1
Size	10 ml
Storage	Stored at 2 - 8°C. Do not freeze. Protect from light.

Advantage

Highly accurate and monodisperse in size, shape

Batch to batch uniformity

Extremely stable

High surface reactivity

High electromagnetic enhancements for SERS and other Raman

Contact Us

Beijing Biotyscience Co. Ltd.

QQ: 499854788

3494243873

WeChat: 13681256816; 15511114213

Email: info@biotyscience.com

Tel: 400-669-8850

15511114213; 13681256816

北京百欧泰生物科技有限公司

Tel: 400-669-8850 Email: info@biotyscience.com

Address: 北京市房山区良乡凯旋大街建设路 18 号