


Gold Nanostars

DESCRIPTION

Gold nanostars have multiple arms with sharp corners, which can generate high electric field amplification at the sharp corners.



Beijing Biotyscience Co. Ltd provides high-quality Gold nanostars. Every batch is completely characterized including size, monodispersity, aggregation, residual chemicals, and concentration.

Storage


Store product away from direct sunlight at 2-8 ° C.

Do NOT freeze. Freezing causes irreversible aggregation of the gold nanoparticles.

When stored as specified the product is stable for at least six months.

PRODUCT INFORMATION

Type	Gold Nanostar
Diameter	70 nm
Surface	CTAB
OD	1
Size	10 ml
Storage	Stored at 2 - 8°C. Do not freeze. Protect from light.



Biotyscience Gold Nanostar

Applications

Biosensors

Surface Plasmon Resonance Enhanced Spectroscopy

Advantage

Exact Sizes

Highly monodisperse sizes and shapes

Long shelf lives

北京百欧泰生物科技有限公司

Tel: 400-669-8850 Email: info@biotyscience.com

Address: 北京市房山区良乡凯旋大街建设路 18 号

High surface reactivity

High electromagnetic enhancements for SERS and other Raman

Contact Us

Beijing Biotyscience Co. Ltd.

QQ: 499854788

3494243873

WeChat: 13681256816; 15511114213

Email: info@biotyscience.com

Tel: 400-669-8850

15511114213; 13681256816