

Palladium Nanoparticles

DESCRIPTION

Biotyscience provides Palladium Nanoparticles of various sizes to meet the needs of different users for different purposes. Biotyscience's Palladium Nanoparticles products use deionized water as the solvent.

The product has a high degree of stability in water and can be evenly dispersed for several months without sedimentation. Each batch of products is synthesized through precise steps and characterized by independent visible/infrared extinction spectroscopy. Different batches of products have excellent repeatability.

PRODUCT INFORMATION

Type	Palladium Nanoparticles
Diameter	3 nm, 20-30 nm
Surface	Citrate (or others)
Concentration	0.05 mg/ml
Size	10 ml
Buffer	Supplied in DI water.
Storage	Stored at 2 - 8°C. Do not freeze. Protect from light.

Storage

Store product away from direct sunlight at 2-8 ° C.

Do NOT freeze. Freezing causes irreversible aggregation of the gold nanoparticles.

When stored as specified the product is stable for at least six months.

Contact Us

Beijing Biotyscience Co. Ltd.

QQ: 499854788

3494243873

WeChat: 13681256816; 15511114213

Email: info@biotyscience.com

Tel: 400-669-8850

15511114213; 13681256816